

For Research Use Only

UFD1L Recombinant antibody

Catalog Number: 83669-1-RR



Basic Information

Catalog Number: 83669-1-RR	GenBank Accession Number: BC005087	Purification Method: Protein A purification
Size: 100ul , Concentration: 1000 µg/ml by Nanodrop;	GeneID (NCBI): 7353	CloneNo.: 230411A2
Source: Rabbit	UNIPROT ID: Q92890	
Isotype: IgG	Full Name: ubiquitin fusion degradation 1 like (yeast)	
Immunogen Catalog Number: AG0954	Calculated MW: 35 kDa	

Applications

Tested Applications:
FC (Intra), ELISA

Species Specificity:
human

Storage

Storage:
Store at -20°C. Stable for one year after shipment.

Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

***** 20ul sizes contain 0.1% BSA**


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)


E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Selected Validation Data



1x10⁶ HEK-293 cells were intracellularly stained with 0.25 ug Ufd1l Recombinant Antibody (83669-1-RR, Clone:230411A2) and APC-Conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)(red), or 0.25 ug Isotype Control (blue). Cells were fixed and permeabilized with Transcription Factor Staining Buffer Kit (PF00011).



1x10⁶ HeLa cells were intracellularly stained with 0.25 ug Ufd1l Recombinant Antibody (83669-1-RR, Clone:230411A2) and APC-Conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)(red), or 0.25 ug Isotype Control (blue). Cells were fixed and permeabilized with Transcription Factor Staining Buffer Kit (PF00011).



Biolayer interferometry (BLI) kinetic assays of 83669-1-RR against Human UFD1L were performed. The affinity constant is below 1 pM.