

For Research Use Only

PAX7 Recombinant antibody, PBS Only (Capture)



Catalog Number: 83679-1-PBS

Basic Information

Catalog Number: 83679-1-PBS	GenBank Accession Number: BC121165	Purification Method: Protein A purification
Size: 100ug , Concentration: 1 mg/ml by Nanodrop;	GeneID (NCBI): 5081	CloneNo.: 240712E6
Source: Rabbit	UNIPROT ID: P23759	
Isotype: IgG	Full Name: paired box 7	
Immunogen Catalog Number: AG25366	Calculated MW: 520 aa, 57 kDa	

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Product Information

83679-1-PBS targets PAX7 as part of a matched antibody pair:

MP00660-1: 83679-1-PBS capture and 83679-3-PBS detection (validated in Cytometric bead array)

MP00660-2: 83679-1-PBS capture and 83679-2-PBS detection (validated in Cytometric bead array)

Unconjugated rabbit recombinant monoclonal antibody in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

This conjugation ready format makes antibodies ideal for use in many applications including: ELISAs, multiplex assays requiring matched pairs, mass cytometry, and multiplex imaging applications. Antibody use should be optimized by the end user for each application and assay.

Storage

Storage:
Store at -80°C.

Storage Buffer:
PBS Only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)


E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00660-1, PAX7 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83679-1-PBS. Detection antibody: 83679-3-PBS. Standard: Ag25366. Range: 0.625-80 ng/mL



Cytometric bead array standard curve of MP00660-2, PAX7 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83679-1-PBS. Detection antibody: 83679-2-PBS. Standard: Ag25366. Range: 2.5-80 ng/mL