

For Research Use Only

TFF2 Recombinant antibody, PBS Only



Catalog Number: 83717-4-PBS

Basic Information

Catalog Number: 83717-4-PBS	GenBank Accession Number: BC032820	Purification Method: Protein A purification
Size: 100ug , Concentration: 1 mg/ml by Nanodrop;	GeneID (NCBI): 7032	CloneNo.: 240661G12
Source: Rabbit	UNIPROT ID: Q03403	
Isotype: IgG	Full Name: trefoil factor 2	
Immunogen Catalog Number: AG4537	Calculated MW: 129 aa, 14 kDa	

Applications

Tested Applications:
WB, ELISA

Species Specificity:
human


Storage

Storage:
Store at -80°C.
Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA) E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Human urine sample was subjected to SDS PAGE followed by western blot with 83717-4-RR (TFF2 antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 83717-4-PBS in a different storage buffer formulation.