

For Research Use Only

TNF-alpha Recombinant antibody, PBS Only



Catalog Number: 83721-7-PBS

Basic Information

Catalog Number: 83721-7-PBS	GenBank Accession Number: NM_000594.4	Purification Method: Protein A purification
Size: 100ug, Concentration: 1 mg/ml by Nanodrop;	GeneID (NCBI): 7124	CloneNo.: 240668B4
Source: Rabbit	ENSEMBL Gene ID: ENSG00000232810	
Isotype: IgG	UNIPROT ID: P01375	
	Full Name: tumor necrosis factor (TNF superfamily, member 2)	
	Calculated MW: 26 kDa	

Applications

Tested Applications:
WB, ELISA

Species Specificity:
human

Storage


Storage:
Store at -80°C.

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA) E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Untreated, LPS and Brefeldin A treated THP-1 cells were subjected to SDS PAGE followed by western blot with 83721-7-RR (TNF-alpha antibody) at dilution of 1:5000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 83721-7-PBS in a different storage buffer formulation.