

For Research Use Only

BBS5 Recombinant antibody, PBS Only (Detector)



Catalog Number: 83722-1-PBS

Basic Information

Catalog Number: 83722-1-PBS	GenBank Accession Number: BC044593	Purification Method: Protein A purification
Size: 100ug , Concentration: 1 mg/ml by Nanodrop;	GeneID (NCBI): 129880	CloneNo.: 240674B5
Source: Rabbit	UNIPROT ID: Q8N3I7	
Isotype: IgG	Full Name: Bardet-Biedl syndrome 5	
Immunogen Catalog Number: AG6153	Calculated MW: 39 kDa	

Applications

Tested Applications:
WB, Cytometric bead array, Indirect ELISA

Species Specificity:
human, mouse, rat

Product Information

83722-1-PBS targets BBS5 as part of a matched antibody pair:

MP00703-1: 83722-2-PBS capture and 83722-1-PBS detection (validated in Cytometric bead array)

Unconjugated rabbit recombinant monoclonal antibody in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

This conjugation ready format makes antibodies ideal for use in many applications including: ELISAs, multiplex assays requiring matched pairs, mass cytometry, and multiplex imaging applications. Antibody use should be optimized by the end user for each application and assay.

Storage

Storage:
Store at -80°C.

Storage Buffer:
PBS Only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)


E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 83722-1-RR (BBS5 antibody) at dilution of 1:5000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 83722-1-PBS in a different storage buffer formulation.



Various lysates were subjected to SDS PAGE followed by western blot with 83722-1-RR (BBS5 antibody) at dilution of 1:5000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 83722-1-PBS in a different storage buffer formulation.



Cytometric bead array standard curve of MP00703-1, BBS5 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83722-2-PBS. Detection antibody: 83722-1-PBS. Standard: Ag6153. Range: 2.5-80 ng/mL.