

For Research Use Only

CBR1 Recombinant antibody, PBS Only



Catalog Number: 83736-4-PBS

Basic Information

Catalog Number: 83736-4-PBS	GenBank Accession Number: BC002511	Purification Method: Protein A purification
Size: 100ug , Concentration: 1 mg/ml by Nanodrop;	GeneID (NCBI): 873	CloneNo.: 240707G1
Source: Rabbit	UNIPROT ID: P16152	
Isotype: IgG	Full Name: carbonyl reductase 1	
Immunogen Catalog Number: AG34830	Calculated MW: 30 kDa	

Applications

Tested Applications:
WB, ELISA

Species Specificity:
human


Storage

Storage:
Store at -80°C.
Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA) E: proteintech@ptglab.com W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 83736-4-RR (CBR1 antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 83736-4-PBS in a different storage buffer formulation.