

For Research Use Only

SETDB1 Recombinant antibody, PBS Only



Catalog Number: 83772-5-PBS

Basic Information

Catalog Number: 83772-5-PBS	GenBank Accession Number: BC009362	Purification Method: Protein A purification
Size: 100ug , Concentration: 1 mg/ml by Nanodrop;	GeneID (NCBI): 9869	CloneNo.: 240805D7
Source: Rabbit	UNIPROT ID: Q15047	
Isotype: IgG	Full Name: SET domain, bifurcated 1	
Immunogen Catalog Number: AG1725	Calculated MW: 143 kDa	

Applications

Tested Applications:
WB, ELISA

Species Specificity:
human, mouse

Storage


Storage:
Store at -80°C.

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA) E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 83772-5-RR (SETDB1 antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 83772-5-PBS in a different storage buffer formulation.