

For Research Use Only

Coagulation Factor III/Tissue Factor Recombinant antibody, PBS Only (Detector)

Catalog Number: 83776-2-PBS



Basic Information

Catalog Number: 83776-2-PBS	GenBank Accession Number: NM_001993.5	Purification Method: Protein A purification
Size: 100ug , Concentration: 1 mg/ml by Nanodrop;	GeneID (NCBI): 2152	CloneNo.: 240784E5
Source: Rabbit	UNIPROT ID: P13726-1	
Isotype: IgG	Full Name: coagulation factor III (thromboplastin, tissue factor)	
	Calculated MW: 33 kDa	

Applications

Tested Applications:
Cytometric bead array, Sandwich ELISA, Indirect ELISA,
Sample test

Species Specificity:
human

Product Information

83776-2-PBS targets Coagulation Factor III/Tissue Factor as part of a matched antibody pair:

MP00732-1: 83776-1-PBS capture and 83776-2-PBS detection (validated in Cytometric bead array, Sandwich ELISA)

Unconjugated rabbit recombinant monoclonal antibody in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

This conjugation ready format makes antibodies ideal for use in many applications including: ELISAs, multiplex assays requiring matched pairs, mass cytometry, and multiplex imaging applications. Antibody use should be optimized by the end user for each application and assay.

Storage


Storage:
Store at -80°C.

Storage Buffer:
PBS Only


For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA) E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Selected Validation Data




Sandwich ELISA standard curve of MP00732-1, Human Coagulation Factor III/Tissue Factor Recombinant Matched Antibody Pair - PBS only. 83776-1-PBS was coated to a plate as the capture antibody and incubated with serial dilutions of standard Eg1052. 83776-2-PBS was HRP conjugated as the detection antibody. Range: 7.8-500 pg/mL



The mean Coagulation Factor III/Tissue Factor concentration was determined to be 79.46 ng/mL in U-87 MG cell extract based on a 3.0 mg/mL extract load.



Serum of sixteen individual healthy human donors was measured. The Coagulation Factor III/Tissue Factor concentration of detected samples was determined to be 44.3 pg/mL with a range of 8.4 - 122.7 pg/mL.



Cytometric bead array standard curve of MP00732-1, Coagulation Factor III/Tissue Factor Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83776-1-PBS. Detection antibody: 83776-2-PBS. Standard: Eg1052. Range: 0.625-80 ng/mL