

For Research Use Only

PRM2 Recombinant antibody, PBS Only



Catalog Number: 83814-1-PBS

Basic Information

Catalog Number: 83814-1-PBS	GenBank Accession Number: BC066338	Purification Method: Protein A purification
Size: 100ug , Concentration: 1 mg/ml by Nanodrop;	GeneID (NCBI): 5620	CloneNo.: 240576B4
Source: Rabbit	UNIPROT ID: P04554	
Isotype: IgG	Full Name: protamine 2	
Immunogen Catalog Number: AG5952	Calculated MW: 13 kDa	

Applications

Tested Applications:
WB, ELISA

Species Specificity:
Human


Storage

Storage:
Store at -80°C.
Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA) E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



human testis tissue were subjected to SDS PAGE followed by western blot with 83814-1-RR (PRM2 antibody) at dilution of 1:20000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 83814-1-PBS in a different storage buffer formulation.