

For Research Use Only

PRMT3 Recombinant antibody, PBS Only (Capture)



Catalog Number: 83969-2-PBS

Basic Information

Catalog Number: 83969-2-PBS	GenBank Accession Number: BC037544	Purification Method: Protein A purification
Size: 100ug , Concentration: 1 mg/ml by Nanodrop;	GeneID (NCBI): 10196	CloneNo.: 241092C9
Source: Rabbit	UNIPROT ID: O60678	
Isotype: IgG	Full Name: protein arginine methyltransferase 3	
Immunogen Catalog Number: AG13466	Calculated MW: 531 aa, 60 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Storage

Storage:
Store at -80°C.
Storage Buffer:
PBS Only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00892-3, PRMT3 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83969-2-PBS. Detection antibody: 83969-3-PBS. Standard: Ag13466. Range: 0.781-100 ng/mL.