

For Research Use Only

RAG2 Recombinant antibody, PBS Only (Capture)



Catalog Number: 83971-1-PBS

Basic Information

Catalog Number: 83971-1-PBS	GenBank Accession Number: BC022397	Purification Method: Protein A purification
Size: 100ug , Concentration: 1 mg/ml by Nanodrop;	GeneID (NCBI): 5897	CloneNo.: 240855C7
Source: Rabbit	UNIPROT ID: P55895	
Isotype: IgG	Full Name: recombination activating gene 2	
Immunogen Catalog Number: AG2393	Calculated MW: 527 aa, 59 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Product Information

83971-1-PBS targets RAG2 as part of a matched antibody pair:

MP00894-3: 83971-1-PBS capture and 83971-3-PBS detection (validated in Cytometric bead array)

Unconjugated rabbit recombinant monoclonal antibody in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

This conjugation ready format makes antibodies ideal for use in many applications including: ELISAs, multiplex assays requiring matched pairs, mass cytometry, and multiplex imaging applications. Antibody use should be optimized by the end user for each application and assay.

Storage

Storage:
Store at -80°C.

Storage Buffer:
PBS Only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00894-3, RAG2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83971-1-PBS. Detection antibody: 83971-3-PBS. Standard: Ag2393. Range: 0.78-100 ng/mL