

For Research Use Only

NFKB2 Recombinant antibody, PBS Only (Capture)



Catalog Number: 84022-1-PBS

Basic Information

Catalog Number: 84022-1-PBS	GenBank Accession Number: BC002844	Purification Method: Protein A purification
Size: 100ug, Concentration: 1 mg/ml by Nanodrop;	GeneID (NCBI): 4791	CloneNo.: 241126D5
Source: Rabbit	ENSEMBL Gene ID: ENSG00000077150	
Isotype: IgG	UNIPROT ID: Q00653	
Immunogen Catalog Number: AG0662	Full Name: nuclear factor of kappa light polypeptide gene enhancer in B-cells 2 (p49/p100)	
	Calculated MW: 97 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Product Information

84022-1-PBS targets NFKB2 as part of a matched antibody pair:

MP00950-3: 84022-1-PBS capture and 84022-4-PBS detection (validated in Cytometric bead array)

Unconjugated rabbit recombinant monoclonal antibody in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

This conjugation ready format makes antibodies ideal for use in many applications including: ELISAs, multiplex assays requiring matched pairs, mass cytometry, and multiplex imaging applications. Antibody use should be optimized by the end user for each application and assay.

Storage

Storage:
Store at -80°C.

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data