

For Research Use Only

SLC9A3R2 Recombinant antibody, PBS Only (Capture)

Catalog Number: 84060-4-PBS



Basic Information

Catalog Number: 84060-4-PBS	GenBank Accession Number: BC106001	Purification Method: Protein A purification
Size: 100ug, Concentration: 1 mg/ml by Nanodrop;	GeneID (NCBI): 9351	CloneNo.: 241215F8
Source: Rabbit	UNIPROT ID: Q15599	
Isotype: IgG	Full Name: solute carrier family 9 (sodium/hydrogen exchanger), member 3 regulator 2	
Immunogen Catalog Number: AG13476	Calculated MW: 337 aa, 37 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Product Information

84060-4-PBS targets SLC9A3R2 as part of a matched antibody pair:

MP00976-1: 84060-4-PBS capture and 84060-1-PBS detection (validated in Cytometric bead array)

MP00976-2: 84060-4-PBS capture and 84060-3-PBS detection (validated in Cytometric bead array)

MP00976-3: 84060-4-PBS capture and 84060-2-PBS detection (validated in Cytometric bead array)

Unconjugated rabbit recombinant monoclonal antibody in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

This conjugation ready format makes antibodies ideal for use in many applications including: ELISAs, multiplex assays requiring matched pairs, mass cytometry, and multiplex imaging applications. Antibody use should be optimized by the end user for each application and assay.

Storage

Storage:
Store at -80°C.

Storage Buffer:
PBS Only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)


E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Selected Validation Data



Cytometric bead array standard curve of MP00976-1, SLC9A3R2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84060-4-PBS. Detection antibody: 84060-1-PBS. Standard: Ag13476. Range: 0.313-40 ng/mL



Cytometric bead array standard curve of MP00976-2, SLC9A3R2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84060-4-PBS. Detection antibody: 84060-3-PBS. Standard: Ag13476. Range: 0.313-40 ng/mL



Cytometric bead array standard curve of MP00976-3, SLC9A3R2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84060-4-PBS. Detection antibody: 84060-2-PBS. Standard: Ag13476. Range: 1.25-40 ng/mL