

For Research Use Only

ZC3H13 Recombinant antibody, PBS Only

Catalog Number: 84228-4-PBS



Basic Information

Catalog Number: 84228-4-PBS	GenBank Accession Number: NM_001076788	Purification Method: Protein A purification
Size: 100ug, Concentration: 1 mg/ml by Nanodrop;	GeneID (NCBI): 23091	CloneNo.: 241516G3
Source: Rabbit	UNIPROT ID: Q5T200	
Isotype: IgG	Full Name: zinc finger CCCH-type containing 13	
Immunogen Catalog Number: AG28909	Calculated MW: 197 kDa	
	Observed MW: 196-240 kDa	

Applications

Tested Applications:
WB, Indirect ELISA

Species Specificity:
human, mouse

Storage


Storage:
Store at -80°C.

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA) E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 84228-4-RR (ZC3H13 antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 84228-4-PBS in a different storage buffer formulation.