

For Research Use Only

Transferrin/TF Recombinant antibody, PBS Only (Capture)

Catalog Number: 84386-4-PBS



Basic Information

Catalog Number: 84386-4-PBS	GenBank Accession Number: BC059367	Purification Method: Protein A purification
Size: 100ug, Concentration: 1 mg/ml by Nanodrop;	GeneID (NCBI): 7018	CloneNo.: 241769D3
Source: Rabbit	UNIPROT ID: P02787	
Isotype: IgG	Full Name: transferrin	
	Calculated MW: 698 aa, 77 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Product Information

84386-4-PBS targets Transferrin/TF as part of a matched antibody pair.

MP01258-1: 84386-4-PBS capture and 84386-2-PBS detection (validated in Cytometric bead array)

Unconjugated rabbit recombinant monoclonal antibody in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

This conjugation ready format makes antibodies ideal for use in many applications including: ELISAs, multiplex assays requiring matched pairs, mass cytometry, and multiplex imaging applications. Antibody use should be optimized by the end user for each application and assay.

Storage

Storage:
Store at -80°C.

Storage Buffer:
PBS Only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01258-1, Transferrin/TF Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84386-4-PBS. Detection antibody: 84386-2-PBS. Standard: Eg0114. Range: 0.781-100 ng/mL