

For Research Use Only

# IL-18R1 Recombinant antibody, PBS Only

Catalog Number: 84390-1-PBS



## Basic Information

Catalog Number:

84390-1-PBS

Size:

100ug , Concentration: 1 mg/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

BC093975

GeneID (NCBI):

8809

UNIPROT ID:

Q13478

Full Name:

interleukin 18 receptor 1

Calculated MW:

541 aa, 62 kDa

Purification Method:

Protein A purification

CloneNo.:

241672F7

## Applications

Tested Applications:

Indirect ELISA

Species Specificity:

human

## Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS only, pH7.3


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Biolayer interferometry (BLI) kinetic assays of 84390-1-PBS against Human IL-18R1 were performed. The affinity constant is 4.36 nM.