

For Research Use Only

# SLC28A3 Recombinant antibody, PBS Only (Capture)

Catalog Number: 84393-2-PBS



## Basic Information

**Catalog Number:**

84393-2-PBS

**Size:**

100ug, Concentration: 1 mg/ml by Nanodrop;

**Source:**

Rabbit

**Isotype:**

IgG

**Immunogen Catalog Number:**

AG12758

**GenBank Accession Number:**

BC093821

**GeneID (NCBI):**

64078

**UNIPROT ID:**

Q9HAS3

**Full Name:**

solute carrier family 28 (sodium-coupled nucleoside transporter), member 3

**Calculated MW:**

691 aa, 77 kDa

**Purification Method:**

Protein A purification

**CloneNo.:**

241615H8

## Applications

**Tested Applications:**

Indirect ELISA

**Species Specificity:**

human

## Product Information

84393-2-PBS targets SLC28A3 as part of a matched antibody pair:

MP01264-1: 84393-2-PBS capture and 84393-1-PBS detection (validated in Cytometric bead array)

Unconjugated rabbit recombinant monoclonal antibody in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

This conjugation ready format makes antibodies ideal for use in many applications including: ELISAs, multiplex assays requiring matched pairs, mass cytometry, and multiplex imaging applications. Antibody use should be optimized by the end user for each application and assay.

## Storage

**Storage:**

Store at -80°C.

**Storage Buffer:**

PBS Only

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

---

## Selected Validation Data