

For Research Use Only

HSPA13 Recombinant antibody, PBS Only (Capture/Detector)

Catalog Number: 84426-1-PBS



Basic Information

Catalog Number: 84426-1-PBS	GenBank Accession Number: BC036370	Purification Method: Protein A purification
Size: 100ug , Concentration: 1 mg/ml by Nanodrop;	GeneID (NCBI): 6782	CloneNo.: 241572B3
Source: Rabbit	UNIPROT ID: P48723	
Isotype: IgG	Full Name: heat shock protein 70kDa family, member 13	
Immunogen Catalog Number: AG3349	Calculated MW: 471 aa, 52 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Product Information

84426-1-PBS targets HSPA13 as part of a matched antibody pair:

MP01302-1: 84426-1-PBS capture and 84426-2-PBS detection (validated in Cytometric bead array)

MP01302-3: 84426-3-PBS capture and 84426-1-PBS detection (validated in Cytometric bead array)

Unconjugated rabbit recombinant monoclonal antibody in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

This conjugation ready format makes antibodies ideal for use in many applications including: ELISAs, multiplex assays requiring matched pairs, mass cytometry, and multiplex imaging applications. Antibody use should be optimized by the end user for each application and assay.

Storage

Storage:
Store at -80°C.

Storage Buffer:
PBS Only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)


E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01302-1, HSPA13 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84426-1-PBS. Detection antibody: 84426-2-PBS. Standard: Ag3349. Range: 0.781-100 ng/mL



Cytometric bead array standard curve of MP01302-3, HSPA13 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84426-3-PBS. Detection antibody: 84426-1-PBS. Standard: Ag3349. Range: 0.781-100 ng/mL