

For Research Use Only

Mouse GM-CSF R alpha Recombinant antibody, PBS Only (Detector)

Catalog Number: 84483-3-PBS



Basic Information

Catalog Number:

84483-3-PBS

Size:

100ug, Concentration: 1 mg/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM_009970.2

GeneID (NCBI):

12982

UNIPROT ID:

Q00941

Full Name:

colony stimulating factor 2 receptor, alpha, low-affinity (granulocyte-macrophage)

Calculated MW:

42kDa

Purification Method:

Protein A purification

CloneNo.:

241889D7

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

mouse

Product Information

84483-3-PBS targets GM-CSF R alpha as part of a matched antibody pair:

MP01319-1: 84483-2-PBS capture and 84483-3-PBS detection (validated in Cytometric bead array)

MP01319-3: 84483-1-PBS capture and 84483-3-PBS detection (validated in Cytometric bead array)

Unconjugated rabbit recombinant monoclonal antibody in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

This conjugation ready format makes antibodies ideal for use in many applications including: ELISAs, multiplex assays requiring matched pairs, mass cytometry, and multiplex imaging applications. Antibody use should be optimized by the end user for each application and assay.

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS Only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)


E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01319-1, MOUSE CSF2RA/CD116 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84483-2-PBS. Detection antibody: 84483-3-PBS. Standard: Eg2016. Range: 0.469-60 ng/mL



Cytometric bead array standard curve of MP01319-3, MOUSE GM-CSF R alpha Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84483-1-PBS. Detection antibody: 84483-3-PBS. Standard: Eg2016. Range: 0.469-60 ng/mL