

For Research Use Only

# CD68 Recombinant antibody, PBS Only (Capture/Detector)

Catalog Number: 84596-3-PBS



## Basic Information

<b>Catalog Number:</b> 84596-3-PBS	<b>GenBank Accession Number:</b> BC015557	<b>Purification Method:</b> Protein A purification
<b>Size:</b> 100ug, Concentration: 1 mg/ml by Nanodrop;	<b>GeneID (NCBI):</b> 968	<b>CloneNo.:</b> 241810G1
<b>Source:</b> Rabbit	<b>UNIPROT ID:</b> P34810	
<b>Isotype:</b> IgG	<b>Full Name:</b> CD68 molecule	
	<b>Calculated MW:</b> 37 kDa	

## Applications

**Tested Applications:**  
Cytometric bead array, Sandwich ELISA, Indirect ELISA,  
Sample test

**Species Specificity:**  
human

## Product Information

84596-3-PBS targets CD68 as part of a matched antibody pair:

MP01429-1: 84596-1-PBS capture and 84596-3-PBS detection (validated in Cytometric bead array)

MP01429-3: 84596-3-PBS capture and 84596-5-PBS detection (validated in Sandwich ELISA)

Unconjugated rabbit recombinant monoclonal antibody in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

This conjugation ready format makes antibodies ideal for use in many applications including: ELISAs, multiplex assays requiring matched pairs, mass cytometry, and multiplex imaging applications. Antibody use should be optimized by the end user for each application and assay.

## Storage

**Storage:**  
Store at -80°C.  
**Storage Buffer:**  
PBS Only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free  
in USA), or 1(312) 455-8498 (outside USA)


E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**


## Selected Validation Data



Sandwich ELISA standard curve of MP01429-3, Human CD68 Recombinant Matched Antibody Pair - PBS only. 84596-3-PBS was coated to a plate as the capture antibody and incubated with serial dilutions of standard Eg2061. 84596-5-PBS was HRP conjugated as the detection antibody. Range: 0.156-10 ng/mL



The mean CD68 concentration was determined to be 10.53 ng/mL in THP-1 cell extract based on a 4.40 mg/mL extract load.



Cytometric bead array standard curve of MP01429-1, CD68 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84596-1-PBS. Detection antibody: 84596-3-PBS. Standard: Eg2061. Range: 0.781-100 ng/mL