

For Research Use Only

STX6 Recombinant antibody, PBS Only (Capture)

Catalog Number: 84627-1-PBS



Basic Information

| | | |
|--|--|---|
| Catalog Number: 84627-1-PBS | GenBank Accession Number: BC009944 | Purification Method: Protein A purification |
| Size: 100ug, Concentration: 1 mg/ml by Nanodrop; | GeneID (NCBI): 10228 | CloneNo.: 242100D8 |
| Source: Rabbit | UNIPROT ID: O43752 | |
| Isotype: IgG | Full Name: syntaxin 6 | |
| Immunogen Catalog Number: AG1284 | Calculated MW: 29 kDa | |

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Product Information

84627-1-PBS targets STX6 as part of a matched antibody pair.

MP01437-1: 84627-1-PBS capture and 84627-2-PBS detection (validated in Cytometric bead array)

Unconjugated rabbit recombinant monoclonal antibody in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

This conjugation ready format makes antibodies ideal for use in many applications including: ELISAs, multiplex assays requiring matched pairs, mass cytometry, and multiplex imaging applications. Antibody use should be optimized by the end user for each application and assay.

Storage

Storage:
Store at -80°C.

Storage Buffer:
PBS Only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MPO1437-1, STX6 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84627-1-PBS. Detection antibody: 84627-2-PBS. Standard: Ag1284. Range: 0.781-100 ng/mL.