

For Research Use Only

TRIM11 Recombinant antibody, PBS Only

Catalog Number:84651-3-PBS



Basic Information

Catalog Number: 84651-3-PBS	GenBank Accession Number: BC011629	Purification Method: Protein A purification
Size: 100ug , Concentration: 1 mg/ml by Nanodrop;	GeneID (NCBI): 81559	CloneNo.: 242107F1
Source: Rabbit	UNIPROT ID: Q96F44	
Isotype: IgG	Full Name: tripartite motif-containing 11	
Immunogen Catalog Number: AG1308	Calculated MW: 52 kDa	
	Observed MW: 53 kDa	

Applications

Tested Applications:
WB, Indirect ELISA

Species Specificity:
human

Background Information

TRIM11(Tripartite motif-containing protein 11) is also named as RNF92 and Belongs to the TRIM/RBCC family. It is a member of the tripartite-motif-containing protein family and is known to destabilize humanin, an inhibitor of Alzheimer-like neuronal insults(PMID: 16904669). This protein has 3 isoforms produced by alternative splicing.

Storage


Storage:
Store at -80°C.

Storage Buffer:
PBS Only


For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA) E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 84651-3-RR (TRIM11 antibody) at dilution of 1:5000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 84651-3-PBS in a different storage buffer formulation.



Biolayer interferometry (BLI) kinetic assays of 84651-3-RR against Human TRIM11 were performed. The affinity constant is 0.25 nM.