

For Research Use Only

PSPC1 Recombinant antibody, PBS Only

Catalog Number: 84657-5-PBS



Basic Information

Catalog Number:

84657-5-PBS

Size:

100ug, Concentration: 1 mg/ml by
Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG10136

GenBank Accession Number:

BC014184

GeneID (NCBI):

55269

UNIPROT ID:

Q8WXF1

Full Name:

paraspeckle component 1

Calculated MW:

523 aa, 59 kDa

Purification Method:

Protein A purification

CloneNo.:

242158C1

Applications

Tested Applications:

IF/ICC, Indirect ELISA

Species Specificity:

human

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS Only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)


E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Immunofluorescent analysis of (4% PFA) fixed HepG2 cells using PSPC1 antibody (84657-5-RR, Clone: 242158C1) at dilution of 1:500 and CoraLite®488-Conjugated Goat Anti-Rabbit IgG(H+L) (SA00013-2), CL594-phalloidin (red). This data was developed using the same antibody clone with 84657-5-PBS in a different storage buffer formulation.



Immunofluorescent analysis of (4% PFA) fixed HepG2 cells using PSPC1 antibody (84657-5-RR, Clone: 242158C1) at dilution of 1:1000 and CoraLite®488-Conjugated Goat Anti-Rabbit IgG(H+L) (SA00013-2), CL594-Phalloidin (red). This data was developed using the same antibody clone with 84657-5-PBS in a different storage buffer formulation.