

For Research Use Only

AGGF1 Recombinant antibody, PBS Only

Catalog Number: 84715-4-PBS



Basic Information

Catalog Number: 84715-4-PBS	GenBank Accession Number: BC029382	Purification Method: Protein A purification
Size: 100ug , Concentration: 1 mg/ml by Nanodrop;	GeneID (NCBI): 55109	CloneNo.: 242264F4
Source: Rabbit	UNIPROT ID: Q8N302	
Isotype: IgG	Full Name: angiogenic factor with G patch and FHA domains 1	
Immunogen Catalog Number: AG2497	Calculated MW: 714 aa, 81 kDa	

Applications

Tested Applications:
IF/ICC, Indirect ELISA

Species Specificity:
human

Background Information

The angiogenic factor gene, AGGF1 (or VG5Q), is identified as a candidate susceptibility gene for Klippel-Trenaunay syndrome (KTS) which is a severe congenital disorder characterized by capillary malformations, venous malformations or varicose veins, and hypertrophy of the affected tissues. AGGF1 protein can bind to endothelial cells and promote cell proliferation. AGGF1 shows strong expression in blood vessels and is secreted as vessel formation is initiated. Regulation of AGGF1 by GATA1 may play roles in endothelial cell biology and angiogenesis.

Storage

Storage:
Store at -80°C.
Storage Buffer:
PBS Only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)


E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Immunofluorescent analysis of (-20°C Ethanol) fixed HeLa cells using AGGF1 antibody (84715-4-RR, Clone: 242264F4) at dilution of 1:250 and CoraLite®488-Conjugated Goat Anti-Rabbit IgG(H+L) (SA00013-2), CL594-Phalloidin (red). This data was developed using the same antibody clone with 84715-4-PBS in a different storage buffer formulation.



Biolayer interferometry (BLI) kinetic assays of 84715-4-RR against Human AGGF1 were performed. The affinity constant is 7.15 nM.