

For Research Use Only

FOSL2 Recombinant antibody, PBS Only

Catalog Number: 84718-5-PBS



Basic Information

Catalog Number: 84718-5-PBS	GenBank Accession Number: BC008899	Purification Method: Protein A purification
Size: 100ug, Concentration: 1 mg/ml by Nanodrop;	GeneID (NCBI): 2355	CloneNo.: 242225C3
Source: Rabbit	UNIPROT ID: P15408	
Isotype: IgG	Full Name: FOS-like antigen 2	
Immunogen Catalog Number: AG8397	Calculated MW: 326 aa, 35 kDa	

Applications

Tested Applications:
IF/ICC, Indirect ELISA

Species Specificity:
human

Storage

Storage:
Store at -80°C.
Storage Buffer:
PBS Only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)


E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Immunofluorescent analysis of (4% PFA) fixed HeLa cells using FOSL2 antibody (84718-5-RR, Clone: 242225C3) at dilution of 1:500 and CoraLite®488-Conjugated Goat Anti-Rabbit IgG(H+L) (SA00013-2), CL594-Phalloidin (red). This data was developed using the same antibody clone with 84718-5-PBS in a different storage buffer formulation.



Biolayer interferometry (BLI) kinetic assays of 84718-5-RR against Human FOSL2 were performed. The affinity constant is 54.4 nM.