

For Research Use Only

# MDM2 Recombinant antibody, PBS Only (Detector)

Catalog Number: 84780-4-PBS



## Basic Information

<b>Catalog Number:</b> 84780-4-PBS	<b>GenBank Accession Number:</b> NM_001145337.2	<b>Purification Method:</b> Protein A purification
<b>Size:</b> 100ug , Concentration: 1 mg/ml by Nanodrop;	<b>GeneID (NCBI):</b> 4193	<b>CloneNo.:</b> 242080C7
<b>Source:</b> Rabbit	<b>UNIPROT ID:</b> Q00987	
<b>Isotype:</b> IgG	<b>Full Name:</b> Mdm2 p53 binding protein homolog (mouse)	
<b>Immunogen Catalog Number:</b> AG27255	<b>Calculated MW:</b> 55 kDa	

## Applications

**Tested Applications:**  
Sandwich ELISA, Indirect ELISA, Sample test

**Species Specificity:**  
human

## Product Information

84780-4-PBS targets MDM2 as part of a matched antibody pair:

MP01554-3: 84780-5-PBS capture and 84780-4-PBS detection (validated in Sandwich ELISA)

Unconjugated rabbit recombinant monoclonal antibody in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

This conjugation ready format makes antibodies ideal for use in many applications including: ELISAs, multiplex assays requiring matched pairs, mass cytometry, and multiplex imaging applications. Antibody use should be optimized by the end user for each application and assay.

## Storage

**Storage:**  
Store at -80°C.

**Storage Buffer:**  
PBS Only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free  
in USA), or 1(312) 455-8498 (outside USA)


E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Sandwich ELISA standard curve of MP01554-3, Human MDM2/HDM2 Recombinant Matched Antibody Pair - PBS only. 84780-5-PBS was coated to a plate as the capture antibody and incubated with serial dilutions of standard Ag27255. 84780-4-PBS was HRP conjugated as the detection antibody. Range: 0.156-10 ng/mL



The mean MDM2/HDM2 concentration was determined to be 12.67 ng/mL in A431 cell extract based on a 1.30 mg/mL extract load.