

For Research Use Only

JMJD5 Recombinant antibody, PBS Only

Catalog Number: 84804-4-PBS



Basic Information

Catalog Number: 84804-4-PBS	GenBank Accession Number: BC027911	Purification Method: Protein A purification
Size: 100ug, Concentration: 1 mg/ml by Nanodrop;	GeneID (NCBI): 79831	CloneNo.: 242097E10
Source: Rabbit	UNIPROT ID: Q8N371	
Isotype: IgG	Full Name: jumonji domain containing 5	
Immunogen Catalog Number: AG15115	Calculated MW: 416 aa, 47 kDa	
	Observed MW: 50-60 kDa	

Applications

Tested Applications:
WB, Indirect ELISA

Species Specificity:
human

Storage


Storage:
Store at -80°C.

Storage Buffer:
PBS Only


For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA) E: proteintech@ptglab.com W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 84804-4-RR (JMJD5 antibody) at dilution of 1:5000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 84804-4-PBS in a different storage buffer formulation.



Biolayer interferometry (BLI) kinetic assays of 84804-4-RR against Human JMJD5 were performed. The affinity constant is 43.3 pM.