

For Research Use Only

UVRAG Recombinant antibody, PBS Only

Catalog Number:84901-1-PBS



Basic Information

Catalog Number: 84901-1-PBS	GenBank Accession Number: BC064837	Purification Method: Protein A purification
Size: 100ug , Concentration: 1 mg/ml by Nanodrop;	GeneID (NCBI): 7405	CloneNo.: 242051C2
Source: Rabbit	UNIPROT ID: Q9P2Y5	
Isotype: IgG	Full Name: UV radiation resistance associated gene	
Immunogen Catalog Number: AG30648	Calculated MW: 699 aa, 78 kDa	
	Observed MW: 90 kDa	

Applications

Tested Applications:
WB, Indirect ELISA

Species Specificity:
human

Background Information

UVRAG is a versatile protein that is involved in the regulation of different cellular pathways implicated in membrane trafficking, and is also involved in the regulation of the COPI-dependent retrograde transport from Golgi and the endoplasmic reticulum by associating with the NRZ complex; Its function is dependent on its binding to phosphatidylinositol 3-phosphate (PtdIns3P) (PMID:16799551; 18552835; 20643123; 24056303; 28306502).

Storage

Storage:
Store at -80°C.

Storage Buffer:
PBS Only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 84901-1-RR (UVRAG antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 84901-1-PBS in a different storage buffer formulation.