

For Research Use Only

SSBP1 Recombinant antibody, PBS Only (Detector)

Catalog Number: 84939-5-PBS



Basic Information

Catalog Number: 84939-5-PBS	GenBank Accession Number: BC000895	Purification Method: Protein A purification
Size: 100ug , Concentration: 1 mg/ml by Nanodrop;	GeneID (NCBI): 6742	CloneNo.: 242519F10
Source: Rabbit	UNIPROT ID: Q04837	
Isotype: IgG	Full Name: single-stranded DNA binding protein 1	
Immunogen Catalog Number: AG2860	Calculated MW: 148 aa, 17 kDa	

Applications

Tested Applications:
Sandwich ELISA, Indirect ELISA, Sample test

Species Specificity:
human

Product Information

84939-5-PBS targets SSBP1 as part of a matched antibody pair:

MP01710-3: 84939-4-PBS capture and 84939-5-PBS detection (validated in Sandwich ELISA)

Unconjugated rabbit recombinant monoclonal antibody in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

This conjugation ready format makes antibodies ideal for use in many applications including: ELISAs, multiplex assays requiring matched pairs, mass cytometry, and multiplex imaging applications. Antibody use should be optimized by the end user for each application and assay.

Storage

Storage:
Store at -80°C.

Storage Buffer:
PBS Only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)


E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Sandwich ELISA standard curve of MP01710-3, Human SSBP1 Recombinant Matched Antibody Pair - PBS only. 84939-4-PBS was coated to a plate as the capture antibody and incubated with serial dilutions of standard Ag2860. 84939-5-PBS was HRP conjugated as the detection antibody. Range: 0.625-40 ng/mL



The mean SSBP1 concentration was determined to be 1,417.48 ng/mL in A549 cell extract based on a 1.50 mg/mL extract load, 3,782.71 ng/mL in A375 cell extract based on a 1.60 mg/mL extract load and 3,185.11 ng/mL in Caco-2 cell extract based on a 2.70 mg/mL extract load.