

For Research Use Only

TFB1M Recombinant antibody, PBS Only

Catalog Number: 84957-4-PBS



Basic Information

Catalog Number:

84957-4-PBS

Size:

100ug, Concentration: 1 mg/ml by
Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG9785

GenBank Accession Number:

BC017788

GeneID (NCBI):

51106

UNIPROT ID:

Q8WVM0

Full Name:

transcription factor B1, mitochondrial

Calculated MW:

346 aa, 40 kDa

Observed MW:

40 kDa

Purification Method:

Protein A purification

CloneNo.:

242524A10

Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

human

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS Only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)


E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 84957-4-RR (TFB1M antibody) at dilution of 1:5000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 84957-4-PBS in a different storage buffer formulation.



Biolayer interferometry (BLI) kinetic assays of 84957-4-RR against Human TFB1M were performed. The affinity constant is 2.99 nM.