

For Research Use Only

# SMARCA5 Recombinant antibody, PBS Only (Capture)

Catalog Number:85019-2-PBS



## Basic Information

**Catalog Number:**

85019-2-PBS

**Size:**

100ug , Concentration: 1 mg/ml by Nanodrop;

**Source:**

Rabbit

**Isotype:**

IgG

**Immunogen Catalog Number:**

AG3759

**GenBank Accession Number:**

BC023144

**GeneID (NCBI):**

8467

**UNIPROT ID:**

O60264

**Full Name:**

SWI/SNF related, matrix associated, actin dependent regulator of chromatin, subfamily a, member 5

**Calculated MW:**

1052 aa, 122 kDa

**Purification Method:**

Protein A purification

**CloneNo.:**

242512A4

## Applications

**Tested Applications:**

Cytometric bead array, Indirect ELISA

**Species Specificity:**

human

## Product Information

85019-2-PBS targets SMARCA5 as part of a matched antibody pair:

MP01767-1: 85019-2-PBS capture and 85019-4-PBS detection (validated in Cytometric bead array)

Unconjugated rabbit recombinant monoclonal antibody in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

This conjugation ready format makes antibodies ideal for use in many applications including: ELISAs, multiplex assays requiring matched pairs, mass cytometry, and multiplex imaging applications. Antibody use should be optimized by the end user for each application and assay.

## Storage

**Storage:**

Store at -80°C.

**Storage Buffer:**

PBS Only

For technical support and original validation data for this product please contact:


T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)

W: [ptglab.com](http://ptglab.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Cytometric bead array standard curve of MP0176-1, SMARCA5 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85019-2-PBS. Detection antibody: 85019-4-PBS. Standard: Ag3759. Range: 0.781-100 ng/mL.