

For Research Use Only

HNRNPR Recombinant antibody, PBS Only (Capture)

Catalog Number:85078-4-PBS



Basic Information

Catalog Number: 85078-4-PBS	GenBank Accession Number: BC001449	Purification Method: Protein A purification
Size: 100ug , Concentration: 1 mg/ml by Nanodrop;	GeneID (NCBI): 10236	CloneNo.: 242734E1
Source: Rabbit	UNIPROT ID: O43390	
Isotype: IgG	Full Name: heterogeneous nuclear ribonucleoprotein R	
Immunogen Catalog Number: AG30339	Calculated MW: 71 kDa	

Applications

Tested Applications:
Indirect ELISA

Species Specificity:
human

Product Information

85078-4-PBS targets HNRNPR as part of a matched antibody pair:

MP01784-3: 85078-4-PBS capture and 85078-2-PBS detection (validated in N/A)

Unconjugated rabbit recombinant monoclonal antibody in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

This conjugation ready format makes antibodies ideal for use in many applications including: ELISAs, multiplex assays requiring matched pairs, mass cytometry, and multiplex imaging applications. Antibody use should be optimized by the end user for each application and assay.

Storage

Storage:
Store at -80°C.

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data