

For Research Use Only

# PRPS2 Recombinant antibody, PBS Only (Capture/Detector)

Catalog Number:85090-1-PBS



## Basic Information

<b>Catalog Number:</b> 85090-1-PBS	<b>GenBank Accession Number:</b> BC119662	<b>Purification Method:</b> Protein A purification
<b>Size:</b> 100ug , Concentration: 1 mg/ml by Nanodrop;	<b>GeneID (NCBI):</b> 5634	<b>CloneNo.:</b> 242657E1
<b>Source:</b> Rabbit	<b>UNIPROT ID:</b> P11908	
<b>Isotype:</b> IgG	<b>Full Name:</b> phosphoribosyl pyrophosphate synthetase 2	
<b>Immunogen Catalog Number:</b> AG25785		

## Applications

**Tested Applications:**  
Cytometric bead array, Indirect ELISA

**Species Specificity:**  
human

## Product Information

85090-1-PBS targets PRPS2 as part of a matched antibody pair:

MP01794-1: 85090-1-PBS capture and 85090-2-PBS detection (validated in Cytometric bead array)

MP01794-2: 85090-3-PBS capture and 85090-1-PBS detection (validated in Cytometric bead array)

Unconjugated rabbit recombinant monoclonal antibody in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

This conjugation ready format makes antibodies ideal for use in many applications including: ELISAs, multiplex assays requiring matched pairs, mass cytometry, and multiplex imaging applications. Antibody use should be optimized by the end user for each application and assay.

## Storage

**Storage:**  
Store at -80°C.  
**Storage Buffer:**  
PBS Only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free  
in USA), or 1(312) 455-8498 (outside USA)


E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Cytometric bead array standard curve of MP01794-1, PRPS2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85090-1-PBS. Detection antibody: 85090-2-PBS. Standard: Ag25785. Range: 0.781-100 ng/mL



Cytometric bead array standard curve of MP01794-2, PRPS2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85090-3-PBS. Detection antibody: 85090-1-PBS. Standard: Ag25785. Range: 0.781-100 ng/mL