

For Research Use Only

# ASCC3 Recombinant antibody, PBS Only (Detector)

Catalog Number: 85130-3-PBS



## Basic Information

<b>Catalog Number:</b> 85130-3-PBS	<b>GenBank Accession Number:</b> BC050681	<b>Purification Method:</b> Protein A purification
<b>Size:</b> 100ug , Concentration: 1 mg/ml by Nanodrop;	<b>GeneID (NCBI):</b> 10973	<b>CloneNo.:</b> 242755C4
<b>Source:</b> Rabbit	<b>UNIPROT ID:</b> Q8N3C0	
<b>Isotype:</b> IgG	<b>Full Name:</b> activating signal cointegrator 1 complex subunit 3	
<b>Immunogen Catalog Number:</b> AG11848	<b>Calculated MW:</b> 111 aa, 13 kDa, 251 kDa	

## Applications

**Tested Applications:**  
Cytometric bead array, Indirect ELISA

**Species Specificity:**  
human

## Product Information

85130-3-PBS targets ASCC3 as part of a matched antibody pair:

MP01828-2: 85130-4-PBS capture and 85130-3-PBS detection (validated in Cytometric bead array)

Unconjugated rabbit recombinant monoclonal antibody in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

This conjugation ready format makes antibodies ideal for use in many applications including: ELISAs, multiplex assays requiring matched pairs, mass cytometry, and multiplex imaging applications. Antibody use should be optimized by the end user for each application and assay.

## Storage

**Storage:**  
Store at -80°C.

**Storage Buffer:**  
PBS only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free  
in USA), or 1(312) 455-8498 (outside USA)

E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Cytometric bead array standard curve of MP01828-2, ASCC3 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85130-4-PBS. Detection antibody: 85130-3-PBS. Standard: Ag11848. Range: 0.781-100 ng/mL.