

For Research Use Only

DNAJC7 Recombinant antibody, PBS Only (Detector)

Catalog Number: 85157-3-PBS



Basic Information

Catalog Number: 85157-3-PBS	GenBank Accession Number: BC011837	Purification Method: Protein A purification
Size: 100ug, Concentration: 1 mg/ml by Nanodrop;	GeneID (NCBI): 7266	CloneNo.: 242875F8
Source: Rabbit	UNIPROT ID: Q99615	
Isotype: IgG	Full Name: DnaJ (Hsp40) homolog, subfamily C, member 7	
Immunogen Catalog Number: AG1572	Calculated MW: 56 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Product Information

85157-3-PBS targets DNAJC7 as part of a matched antibody pair:

MP01884-1: 85157-2-PBS capture and 85157-3-PBS detection (validated in Cytometric bead array)

MP01884-2: 85157-1-PBS capture and 85157-3-PBS detection (validated in Cytometric bead array)

Unconjugated rabbit recombinant monoclonal antibody in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

This conjugation ready format makes antibodies ideal for use in many applications including: ELISAs, multiplex assays requiring matched pairs, mass cytometry, and multiplex imaging applications. Antibody use should be optimized by the end user for each application and assay.

Storage

Storage:
Store at -80°C.

Storage Buffer:
PBS Only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)


E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01884-1, DNJC7 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85157-2-PBS. Detection antibody: 85157-3-PBS. Standard: Ag1572. Range: 0.781-100 ng/mL



Cytometric bead array standard curve of MP01884-2, DNJC7 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85157-1-PBS. Detection antibody: 85157-3-PBS. Standard: Ag1572. Range: 0.781-100 ng/mL