

For Research Use Only

SHC4 Recombinant antibody, PBS Only

Catalog Number:85181-4-PBS



Basic Information

| | | |
|-------------------------------------------------------------|------------------------------------------------------------------------------|-------------------------------------------------------|
| Catalog Number: 85181-4-PBS | GenBank Accession Number: BC033907 | Purification Method: Protein A purification |
| Size: 100ug , Concentration: 1 mg/ml by Nanodrop; | GeneID (NCBI): 399694 | CloneNo.: 242805E3 |
| Source: Rabbit | UNIPROT ID: Q6S5L8 | |
| Isotype: IgG | Full Name: SHC (Src homology 2 domain containing) family, member 4 | |
| Immunogen Catalog Number: AG3336 | Calculated MW: 630 aa, 69 kDa | |
| | Observed MW: 75 kDa | |

Applications

Tested Applications:
WB, Indirect ELISA

Species Specificity:
human

Storage

Storage:
Store at -80°C.

Storage Buffer:
PBS Only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



fetal human brain tissue were subjected to SDS PAGE followed by western blot with 85181-4-RR (SHC4 antibody) at dilution of 1:15000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 85181-4-PBS in a different storage buffer formulation.