

For Research Use Only

ALDH1B1 Recombinant antibody, PBS Only

Catalog Number: 85426-3-PBS



Basic Information

Catalog Number: 85426-3-PBS	GenBank Accession Number: BC001619	Purification Method: Protein A purification
Size: 100ug , Concentration: 1 mg/ml by Nanodrop;	GeneID (NCBI): 219	CloneNo.: 242944F5
Source: Rabbit	UNIPROT ID: P30837	
Isotype: IgG	Full Name: aldehyde dehydrogenase 1 family, member B1	
Immunogen Catalog Number: AG7144	Calculated MW: 57 kDa	

Applications

Tested Applications:
Indirect ELISA

Species Specificity:
human

Product Information

85426-3-PBS targets ALDH1B1 as part of a matched antibody pair:

MP01995-2: 85426-1-PBS capture and 85426-3-PBS detection (validated in Cytometric bead array)

Unconjugated rabbit recombinant monoclonal antibody in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

This conjugation ready format makes antibodies ideal for use in many applications including: ELISAs, multiplex assays requiring matched pairs, mass cytometry, and multiplex imaging applications. Antibody use should be optimized by the end user for each application and assay.

Storage

Storage:
Store at -80°C.

Storage Buffer:
PBS only, pH7.3

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

