

For Research Use Only

# CYP17A1 Recombinant antibody, PBS Only (Capture)

Catalog Number:85690-6-PBS



## Basic Information

<b>Catalog Number:</b> 85690-6-PBS	<b>GenBank Accession Number:</b> BC062997	<b>Purification Method:</b> Protein A purification
<b>Size:</b> 100ug , Concentration: 1 mg/ml by Nanodrop;	<b>GeneID (NCBI):</b> 1586	<b>CloneNo.:</b> 243076G4
<b>Source:</b> Rabbit	<b>UNIPROT ID:</b> P05093	
<b>Isotype:</b> IgG	<b>Full Name:</b> cytochrome P450, family 17, subfamily A, polypeptide 1	
<b>Immunogen Catalog Number:</b> AG5899	<b>Calculated MW:</b> 57 kDa	

## Applications

**Tested Applications:**  
Sandwich ELISA, Indirect ELISA, Sample test

**Species Specificity:**  
human

## Product Information

85690-6-PBS targets CYP17A1 as part of a matched antibody pair:

MP02011-3: 85690-6-PBS capture and 85690-2-PBS detection (validated in Sandwich ELISA)

Unconjugated rabbit recombinant monoclonal antibody in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

This conjugation ready format makes antibodies ideal for use in many applications including: ELISAs, multiplex assays requiring matched pairs, mass cytometry, and multiplex imaging applications. Antibody use should be optimized by the end user for each application and assay.

## Storage

**Storage:**  
Store at -80°C.  
**Storage Buffer:**  
PBS Only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)


E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Sandwich ELISA standard curve of MP02011-3, Human CYP17A1 Recombinant Matched Antibody Pair - PBS only. 85690-6-PBS was coated to a plate as the capture antibody and incubated with serial dilutions of standard Ag5899. 85690-2-PBS was HRP conjugated as the detection antibody. Range: 1.56-100 ng/mL



The mean CYP17A1 concentration was determined to be 93.13 ng/mL in A549 cell extract based on a 3.10 mg/mL extract load and 105.10 ng/mL in HepG2 cell extract based on a 1.80 mg/mL extract load.