

For Research Use Only

Nectin 3 Recombinant antibody

Catalog Number: 85833-4-RR



Basic Information

Catalog Number:

85833-4-RR

Size:

100ul, Concentration: 1000 µg/ml by

Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

EG2657

GenBank Accession Number:

NM_015480.3

GeneID (NCBI):

25945

UNIPROT ID:

Q9NQS3-1

Full Name:

poliovirus receptor-related 3

Calculated MW:

61 kDa

Observed MW:

70-80 kDa

Purification Method:

Protein A purification

CloneNo.:

242804D5

Recommended Dilutions:

WB 1:500-1:2000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human, mouse

Positive Controls:

WB : HCT 116 cells, mouse testis tissue

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol, pH7.3

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 85833-4-RR (PVRL3 antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours.