

For Research Use Only

ACADVL Recombinant antibody, PBS Only

Catalog Number: 85889-1-PBS



Basic Information

Catalog Number:

85889-1-PBS

Size:

100ug , Concentration: 1 mg/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG6016

GenBank Accession Number:

BC012912

GeneID (NCBI):

37

UNIPROT ID:

P49748

Full Name:

acyl-Coenzyme A dehydrogenase, very long chain

Calculated MW:

70 kDa

Observed MW:

60 kDa

Purification Method:

Protein A purification

CloneNo.:

250199C12

Applications

Tested Applications:

WB, IF/ICC, Indirect ELISA

Species Specificity:

human, mouse

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS only, pH7.3


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)


E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 85889-1-RR (ACADVL antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 85889-1-PBS in a different storage buffer formulation.



Immunofluorescent analysis of (4% PFA) fixed A549 cells using ACADVL antibody (85889-1-RR, Clone: 250199C12) at dilution of 1:500 and CoraLite® 488-Conjugated Goat Anti-Rabbit IgG(H+L) (SA00013-2), CL594-Phalloidin (red). This data was developed using the same antibody clone with 85889-1-PBS in a different storage buffer formulation.