

For Research Use Only

FOXJ1 Recombinant antibody, PBS Only

Catalog Number: 86070-1-PBS



Basic Information

Catalog Number:

86070-1-PBS

Size:

100ug, Concentration: 1 mg/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG28294

GenBank Accession Number:

BC046460

GeneID (NCBI):

2302

UNIPROT ID:

Q92949

Full Name:

forkhead box J1

Calculated MW:

45 kDa

Observed MW:

50 kDa

Purification Method:

Protein A purification

CloneNo.:

250724E9

Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

human

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS only, pH7.3


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 86070-1-RR (FOXJ1 antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 86070-1-PBS in a different storage buffer formulation.