

For Research Use Only

# Anti-Human Fas/CD95 Rabbit Recombinant Antibody

Catalog Number: 98034-2-RR



## Basic Information

<b>Catalog Number:</b> 98034-2-RR	<b>GenBank Accession Number:</b> NM_000043.6	<b>Purification Method:</b> Protein A purification
<b>Size:</b> 100ug, 1 mg/ml	<b>GeneID (NCBI):</b> 355	<b>CloneNo.:</b> 240340G11
<b>Source:</b> Rabbit	<b>UNIPROT ID:</b> P25445-1	
<b>Isotype:</b> IgG	<b>Full Name:</b> Fas (TNF receptor superfamily, member 6)	
	<b>Calculated MW:</b> 38 kDa	

## Applications

**Tested Applications:**  
FC

**Species Specificity:**  
human

## Background Information

FAS, also named as CD95, APO-1, APT1, FAS1 and TNFRSF6, is a receptor for TNFSF6/FASLG. It is a cell surface receptor belonging to the TNF receptor superfamily, can mediate apoptosis by ligation with an agonistic anti-Fas antibody or Fas ligand. Stimulation of Fas results in the aggregation of its intracellular death domains, leading to the formation of the death-inducing signaling complex (DISC). FAS-mediated apoptosis may have a role in the induction of peripheral tolerance, in the antigen-stimulated suicide of mature T-cells, or both.

## Storage

**Storage:**  
Store at 2-8°C. Stable for one year after shipment.

**Storage Buffer:**  
PBS with 0.09% sodium azide.


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)


E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**


## Selected Validation Data



1x10<sup>6</sup> human PBMCs were surface stained with 0.25 ug Anti-Human Fas/CD95 Rabbit Recombinant Antibody (98034-2-RR, Clone: 240340G11) or Isotype Control and PE-Conjugated AffiniPure Goat Anti-Rabbit IgG(H+L). Cells were co-stained with CoraLite® Plus 647 Anti-Human CD45RO (UCHL1) (CL647-65150, Clone: UCHL1). Cells were not fixed.



1x10<sup>6</sup> human PBMCs were surface stained with 0.25 ug Anti-Human Fas/CD95 Rabbit Recombinant Antibody (98034-2-RR, Clone: 240340G11) (red) or Isotype Control (blue), and PE-Conjugated AffiniPure Goat Anti-Rabbit IgG(H+L). Cells were not fixed.



Bi-layer interferometry (BLI) kinetic assays of 98034-2-RR against Human Fas/CD95 were performed. The affinity constant is 0.191 nM.