

For Research Use Only

Anti-Human VE-cadherin/CD144 Rabbit Recombinant Antibody

Catalog Number: 98071-1-RR



Basic Information

Catalog Number: 98071-1-RR	GenBank Accession Number: NM_001795.5	Purification Method: Protein A purification
Size: 100ug , 1000 µg/ml	GeneID (NCBI): 1003	CloneNo.: 240755B2
Source: Rabbit	UNIPROT ID: P33151-1	
Isotype: IgG	Full Name: cadherin 5, type 2 (vascular endothelium)	
	Calculated MW: 88 kDa	

Applications

Tested Applications:
FC

Species Specificity:
human

Background Information

Cadherins are a family of transmembrane glycoproteins that mediate calcium-dependent cell-cell adhesion and play an important role in the maintenance of normal tissue architecture. Vascular endothelial cadherin (VE-cadherin), also known as Cadherin-5 (CDH5) or CD144, is a member of the type II classical cadherin family of cell adhesion proteins (PMID: 21269602). VE-cadherin is expressed specifically in endothelial cells and mediates homophilic adhesion in the vascular endothelium (PMID: 1522121; 8555485; 21269602). VE-cadherin plays a role in the organization of lateral endothelial junctions and in the control of permeability properties of vascular endothelium (PMID: 1522121). VE-cadherin has also been shown to be required for angiogenesis (PMID: 16473763; 18162609). The calculated molecular weight of VE-cadherin is 88 kDa and the apparent molecular weight of 120-140 kDa is higher due to post-translational glycosylation and phosphorylation (PMID: 10460833; 29894844). Full-length VE-cadherin can be proteolytically cleaved to generate a fragment of 90-100 kDa (PMID: 9786462; 22064597).

Storage


Storage:
Store at 2 - 8°C. Stable for one year after shipment.

Storage Buffer:
PBS with 0.09% sodium azide, pH 7.3.


For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA) E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



1×10^6 HUVEC cells were stained with 0.25 μ g Anti-Human VE-cadherin/CD144 Rabbit Recombinant Antibody (98071-1-RR, Clone:240755B2) (red), or 0.25 μ g Isotype Control (blue). Cells were not fixed.



Biolayer interferometry (BLI) kinetic assays of 98071-1-RR against Human VE-cadherin/CD144 were performed. The affinity constant is below 1 pM.