

For Research Use Only

rat anti Myc-tag Monoclonal antibody (9E1)



www.ptglab.com

Catalog Number: 9e1

38 Publications

Basic Information

Catalog Number:
9e1

Clone No.:
9E1

Applications:
WB, IF, ELISA, IP

Host/IsoType:
Rat/IgG1

Conjugate:
Unconjugated

Type:
Primary Antibody

Class:
Monoclonal

RRID:
AB_2631398

Purification Method:
Protein G purified

Description

Rat monoclonal [9E1] to Myc-tag

Specificity/Target

Myc-tag (EQKLISEEDL) at the N-terminus, C-terminus, or internal site of the fusion protein.

Recommended Dilution

Western Blot: 1:1,000
Immunofluorescence: 1:500
ELISA: 1:1,000
Immunoprecipitation: 5 µg/ IP

Storage

Storage:
Shipped on blue ice. Store at +4°C. Stable for one year.


Storage Buffer:
PBS, preservative: 0.09% sodium azide

For technical support and original validation data for this product please contact


T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or E: proteintech@ptglab.com
1(312) 455-8498 (outside USA) W: www.ptglab.com

This product is exclusively available under
Proteintech Group brand and is not available
to purchase from any other manufacturer.

Selected Validation Data



Western Blot analysis of cell extract from HEK293T cells expressing PCNA-myc



Immunofluorescence of HeLa cells transiently expressing MAG-myc