

For Research Use Only

FcZero-rAb™ APC Anti-Human CD8a (OKT8) Rabbit IgG Recombinant Antibody

Catalog Number: APC-FcA65506



Basic Information

Catalog Number: APC-FcA65506	GenBank Accession Number: BC025715	Purification Method: Protein A purification
Size: 100tests, 5 ul/test	GeneID (NCBI): 925	CloneNo.: OKT8
Source: Rabbit	ENSEMBL Gene ID: ENSG00000153563	Excitation/Emission maxima wavelengths: 650 nm / 660 nm
Isotype: IgG	Full Name: CD8a molecule	
	Calculated MW: 235 aa, 26 kDa	

Applications

Tested Applications:
FC

Species Specificity:
human

Background Information

CD8 is a transmembrane glycoprotein composed of two disulfide-linked chains. It can be present as a homodimer of CD8a or as a heterodimer of CD8a and CD8β (PMID: 3264320; 8253791). CD8 is found on most thymocytes. The majority of class I-restricted T cells express mostly the CD8aβ heterodimer while CD8aa homodimers alone have been found on some gut intraepithelial T cells, on some T cell receptor (TCR) γδ T cells and on NK cells (PMID: 2111591; 1831127; 8420975). CD8 acts as a co-receptor that binds to MHC class-I and participates in cytotoxic T-cell activation (PMID: 8499079). During T cell development, CD8 is required for positive selection of CD4-/CD8+ T cells (PMID: 1968084).

Storage


Storage:
Store at 2-8°C. Avoid exposure to light. Stable for one year after shipment.

Storage Buffer:
PBS with 0.09% sodium azide and 0.5% BSA.

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA) E: proteintech@ptglab.com W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



1x10⁶ human PBMCs were surface stained with APC Anti-Human CD4 and 5 ul PE Anti-Human CD8a (OKT8) Rabbit IgG Recombinant Antibody (PE-FcA65506, Clone: OKT8) or PE Rabbit IgG Isotype Control Recombinant Antibody (PE-FcA98136, Clone: 240953C9). Cells were not fixed.