

For Research Use Only

FcZero-rAb™ APC Anti-Human ICAM-1/CD54 Rabbit Recombinant Antibody

Catalog Number: APC-FcA98011



Basic Information

Catalog Number: APC-FcA98011	GenBank Accession Number: BC015969	Purification Method: Protein A purification
Size: 100tests , 5 ul/test	GeneID (NCBI): 3383	CloneNo.: 240646A11
Source: Rabbit	ENSEMBL Gene ID: ENSG00000090339	Excitation/Emission maxima wavelengths: 650 nm / 660 nm
Isotype: IgG	Full Name: intercellular adhesion molecule 1	
	Calculated MW: 90 kDa	

Applications

Tested Applications:
FC

Species Specificity:
human

Background Information

ICAM-1 (CD54) is a 90-kDa transmembrane glycoprotein of the immunoglobulin superfamily and is critical for the firm attachment and transmigration of leukocytes out of blood vessels and into tissues (PMID: 19307690). ICAM-1 is expressed by several cell types, typically on endothelial cells and cells of the immune system, and its expression can be up-regulated by various stimuli, including TNF- α , INF- γ , IL-1 and thrombin (PMID: 3086451; 9694714; 15979056). It is a ligand for LFA-1 and Mac-1, serves as a receptor for rhinovirus, and is one of several receptors used by Plasmodium falciparum (PMID: 2566624; 2538244; 2475784).

Storage

Storage:
Store at 2-8°C. Avoid exposure to light. Stable for one year after shipment.

Storage Buffer:
PBS with 0.09% sodium azide and 0.5% BSA.


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



1x10⁶ human PBMCs were surface stained with PE Anti-Human CD19 and 5 ul APC Anti-Human ICAM-1/CD54 Rabbit RecAb (APC-FcA98011, Clone:240646A11), or APC Rabbit IgG Isotype Control RecAb (APC-FcA98136, Clone: 240953C9). Cells were not fixed.