

For Research Use Only

FcZero-rAb™ APC Anti-Human CEACAM5/8 (CD66b/e) Rabbit Recombinant Antibody

Catalog Number: APC-FcA98012



Basic Information

Catalog Number: APC-FcA98012	GenBank Accession Number: BC026263	Purification Method: Protein A purification
Size: 100tests , 5 ul/test	GeneID (NCBI): 1088	CloneNo.: 230353B7
Source: Rabbit	Full Name: carcinoembryonic antigen-related cell adhesion molecule 8	Excitation/Emission maxima wavelengths: 650 nm / 660 nm
Isotype: IgG	Calculated MW: 38 kDa	
Immunogen Catalog Number: EG32133		

Applications

Tested Applications:
FC

Species Specificity:
human

Background Information

Carcinoembryonic antigen-related cell adhesion molecules (CEACAMs) belong to the immunoglobulin superfamily of cell adhesion molecules. CEACAMs play major roles in cell-cell and cell-ECM adhesion and have been implicated in controlling cell proliferation, angiogenesis, and tissue remodeling (PMID: 23021083). Carcinoembryonic antigen (CEA), also known as CEACAM5 or CD66e mainly serves as a cell adhesion molecule mediating intercellular contact by both homophilic and heterophilic binding (PMID: 21731662). CEACAM8, also known as CD66b, CGM6 or NCA-95, is a glycosylphosphatidylinositol (GPI)-linked glycoprotein expressed on granulocytes and is involved in regulating adhesion and activation of eosinophils (PMID: 8104998; 18056392; 18056392). This antibody is raised against 35-319 aa of human CEACAM8/CD66b and can cross-react with CEACAM5/CD66e.

Storage


Storage:
Store at 2-8°C. Avoid exposure to light. Stable for one year after shipment.

Storage Buffer:
PBS with 0.09% sodium azide and 0.5% BSA, pH7.3

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



1x10⁶ human peripheral blood leukocytes were surface stained with 5 ul APC Anti-Human CEACAM5/8 (CD66b/e) Rabbit RecAb (APC-FcA98012, Clone: 230353B7), or 5 ul APC Rabbit IgG Isotype Control RecAb (APC-FcA98136, Clone: 240953C9). Cells were not fixed.