

FOR IN VITRO RESEARCH USE ONLY.  
NOT FOR USE IN HUMANS OR ANIMALS.

# Recombinant human ACE2 protein



## Basic Information

**Catalog Number:**

Ag15455

**Size:**

50 µg

**Form:**

Available lyophilized

**Species:**

human

**Expression Source:**

*e coli*-derived, PGEX-4T, with N-terminal GST.

**Biological Activity:**

Not tested

**Endotoxin Level:**

Please contact the lab for more information

**Peptide Sequence:**

LRNGANEGFHEAVGEIMLSAATPKHLKSI GLLSPDFQ  
EDNETEINFLKQALTIVGTLPFITYMLEKWRWMVFKG  
EIPKDQWMKKWWEMKREIVGVVEVPVPHDETYCDPA  
SLFHVSNDSYFIRYYTRTLYQFQFQEALCQAAKHEGP  
LHKCDISNSTEAGQKLFNMLRLGKSEPWTALENVVG  
AKNMNVRPLLNYFEPLFTWLKDQNKNSFVGWSTDWS  
PYADQSIKVRISLKSALGDKAYEWNDEMFLFRSSVA  
YAMRQYFLKVKNQMLFG EEDVRVANLKPRI SFNFFV  
TAPKNVSDIIPRTEVEKAI RMSRSRINDA FRLNDNSLEF  
LGIQPTLGPNNQPPVSIWLI

(392-744 aa encoded by BC048094)

## Selected Validation Data

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522)  
(toll free in USA), or 1(312) 455-8498  
(outside USA)

E : proteintech@ptglab.com W:  
ptglab.com

This product is exclusively available under Proteintech  
Group brand and is not available to purchase from any  
other manufacturer.