

FOR IN VITRO RESEARCH USE ONLY.  
NOT FOR USE IN HUMANS OR ANIMALS.

# Recombinant human PPP1R3A protein



## Basic Information

### Catalog Number:

Ag18013

### Size:

50 µg

### Form:

Available lyophilized

### Species:

human

### Expression Source:

*e coli*-derived, PGEX-4T, with N-terminal GST.

### Biological Activity:

Not tested

### Endotoxin Level:

Please contact the lab for more information

### Peptide Sequence:

FDPHEGRNDDSHYTLQQRDTVGVIYDNDFEKESRLGI  
CNVRVDEMEKEETMSMYNPRKTHDREKCGTGNITSV  
EESWWITEYQKATSKLDLQLGMLPTDKTVFSENDRH  
RQVQELSKKTDSDAIVHSAFNSDTNRAPQNSPFSK  
HHTEISVSTNEQAIAVENAVTTMASQPISTKSENICN  
STREIQGIEKHPYPESKPEEVSRSSTGIVTSGSRKERCIG  
QIFQTEEYSVEKSLGPMILINKPLENMEEARHENEGLV  
SSGQSLYTSGEKESDSSASTSLPVEESQAQGNESLFS  
KYTNISKIPYFLFLFLITVYHYDLMIGLTFVLSLSWL  
WEEGRQKESVKKK

(773-1122 aa encoded by BC126453)

## Selected Validation Data

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522)  
(toll free in USA), or 1(312) 455-8498  
(outside USA)

E : proteintech@ptglab.com W:  
ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.