

FOR IN VITRO RESEARCH USE ONLY.
NOT FOR USE IN HUMANS OR ANIMALS.

RECOMBINANT HUMAN IRF7



Basic Information

Catalog Number:

Ag18059

Size:

50 µg

Form:

Available lyophilized

Species:

human

Expression Source:

E. coli-derived, PGEX-4T, with N-terminal GST.

Biological Activity:

Not tested

Endotoxin Level:

Please contact the lab for more information

Validated Application:

Blocking peptide

Peptide Sequence:

PGPFLAHTHAGLQAPGPLPAPAGDKGDLQLQAVQQS
CLADHLLTASWGADPVPTKAPGEGQEGLPLTGACAG
GPGLPAGELYGWAVETTPSPGPQPAALTTGEAAPE
SPHQAEPYLSPSPSACTAVQEPSGALDVTIMYKGR
VLQKVVGHPSCTFLYGPPDPAVRATDPQQVAFPS
ELPDQKQLRYTEELLRHVAPGLHLELRGQLWARRM
KCKVYWEVGGPPGSASPSTPACLLPRNCDTPIFDFR
FFQELVEFRARQRRGSPRYIYLGFGQDLSAGRPKE
LVLVKLEPWLCRVHLEGTQREGVSSLDSSSLCLSSA
NSLYDDIECFLMELEQPA

(168-516 aa encoded by BC136555)

Selected Validation Data

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA),
or 1(312) 455-8498 (outside USA)

E : proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.